

RegenLab

Linee Guida in

MEDICINA RIGENERATIVA

ISO 9001-2000
ISO 13485-2003
ISO 14644 (GLP)

SWISSmedic

L'Azienda

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ABSTRACT

SHOULDER SURGEON AND AUTOLOGOUS CELLULAR REGENERATION (ACR)- FROM BENCH TO BED: PART ONE- THE LINK BETWEEN THE HUMAN FIBROBLAST, CONNECTIVE TISSUE DISORDERS AND SHOULDER.

PHD. DPHIL (OXON), FCS., FRCS.

*¹RESEARCH FELLOW, DIVISION OF ANATOMY AND HISTOLOGY,
ACADEMIC DEPARTMENT OF BIOMEDICAL SCIENCES,
FACULTY OF HEALTH SCIENCES,
UNIVERSITY OF STELLENBOSCH, SOUTH AFRICA.*


*¹DIABETES MELLITUS DISCOVERY RESEARCH PLATFORM, MEDICAL RESEARCH COUNCIL, PAROW, SOUTH AFRICA
et al.*

DATE: 29 July 2007(version 1)

The process of autologous cell regeneration (ACR) is a facet of cell therapy and regenerative medicine. It is initiated when activated autologous platelet rich plasma alone (PRP), containing cytokines or growth factors, is injected into the dermis or other structure in order to initiate a regenerative or antiaging process. The recipients resident cells at the target zone are activated by biologically active growth factors, derived from the activated platelets in the PRP through a paracrine effect. The platelet gel that contains the fibrin and clumped platelets, releases growth factors that influences activation of macrophages, and stem cells in the recipient site. Tissue regeneration is facilitated by stem cell proliferation and differentiation. The PRP that has sealant and wound healing properties, may shift the wound healing cascade to the left, thereby speeding up tissue regeneration and remodelling by the use of the recipients own plasma.

Both the fibroblast and myofibroblast play key roles in the wound healing cascade. The fibroblast, of mesenchymal origin, plays a pivotal role in the formation of the extracellular matrix and deposition of collagen. A failure of organised fibroblast function results in important and disabling disease processes and conditions such as chronic ischaemic heart disease and remodelling of the heart, lung fibrosis, fibromatosis, solar aged-face, keloids, hypertrophic scar formation, nodular fasciitis, inguinal hernia, Dupuytren's disease and scleroderma .More recent additions include frozen shoulder and the captured shoulder. In aesthetic medicine, one of the prime functions of fractional photothermolysis with lasers, intense pulsed light and radiofrequency devices, is to stimulate dermal fibroblast proliferation, by thermal energy thereby increasing collagen deposition which enhances facial rejuvenation. Platelet –rich plasma (REGENLAB PRP, REGENLAB-ACR) possesses unique growth factors that stimulate, fibroblast, keratinocyte and myoblasts ex vivo in tissue culture, allowing 3 dimensional cell proliferation within the fibrin gel. PRP can be used as a cell carrier (ie keratinocytes, fibroblasts), may enhance cell retention at the point of treatment. REGENKIT is authorized for human use, ISO and CE marked. In this overview fibroblast morphology, tissue culture and cell biology relevant to the shoulder surgeon is reviewed.

WOUND HEALING

ABSTRACT

CAN WE INCREASE THE RATE OF WOUND HEALING BY USING AUTOLOGOUS PLASMA SOLUTION OBTAINED WITH THE REGEN PRP KIT?

PD Dr med. Wassim Raffoul, et al.

The aim of this clinical research was to evaluate the potential for wound healing of autologous plasma preparation enriched with thrombocytes obtained with the Regen PRP (platelet rich plasma) kit. The Regen PRP kit was evaluated for its simplicity and rapidity of use in real situations such as « point of care » pre- or perioperative procedures related to plastic surgery. The kit allows the preparation of an autologous solution of plasma enriched with thrombocytes, having a minimal amount of erythrocytes, in order to obtain a « white clot » in the wound bed. This PRP preparation, derived from the blood of the patient, can be used in the field of surgery, for treatment of severe burns and for the transplantation of cells and tissues.

It is well known that leukocytes and thrombocytes, imbedded in the fibrin mesh surrounding a wound, synthesise and release all molecules required for the healing process following a specific time sequence. Using a solution of plasma enriched with thrombocytes prepared with Regen PRP kit, an high concentration of thrombocytes is obtained which results in an increase of amount of growth factors and signal molecules available for an optimal and faster healing process. Also, the PRP preparation contains all the supplementary factors (vitamins, essential elements etc.) Because all constituents retrieved in the PRP are from autologous origin, the risk of a local immunological reaction is ruled out.

ABSTRACT

INNOVATION IN DERMATOLOGY : The use of autologous platelet-enriched plasma in three cases

M. Rave¹, C.H. BraccoNolin CH¹, J. Berthou¹, M.H. Choulet¹, A. Grumblat¹, S. Cairey-Remonnay

Pharmacy Department, 2Dermatology Department, Besançon University Hospital, France

BACKGROUND

Emerging cellular therapies such as platelet-rich plasma (PRP) provides wound management options. Local application of a solution of plasma enriched with thrombocytes derived from the patient's blood intends to increase the wound repair rate and reduce pain. The efficacy of the treatment has been proven in surgical applications, in the treatment of severe burns and in transplantation of cells and tissues. Our study deals with leg ulcers. Protocol was performed with ReGen Lab product : ReGenPRP-Kit.

METHODS

The treatment involved one application of autologous PRP to the wound every 15 days. Dressings were changed every 2 days using Jelonet® or Melolin® according to local protocol. The total wound surface area was measured before the first treatment (T0), after two treatments (T2) and after five treatments (T5).

RESULTS

Three patients were treated : 2 for chronic ulcer and 1 for perforated ulcer of the foot (table). In each case, there was rapid epidermisation after the first application of ReGenPRP.

	Age (years)	Sex	Risk factors	Location of lesion and pathology	No. of lesions	No. of treatments	Total surface area of lesion (cm ²)		
							T0	T2	T5
Patient 1	83	Female	Lower limb obstructive arteriopathy Rheumatoid arthritis Knee prosthesis	Mixed chronic ulcers on the two lower limbs, with tendons exposed	5	8	97.6	72.2	20.7
Patient 2	80	Female	Post-phlebitis venous insufficiency Knee and hip prosthesis Allergy	Chronic venous ulcer	2	8	109.6	56.3	42.8
Patient 3	58	Male	Underlying chronic osteitis Diabetes Chronic alcoholism	Perforated ulcers of the foot	1	9	12.3	8.7	7.4

CONCLUSION

This new strategy provides varied results in patients ; some of whom with encouraging results were hospitalised for several months after failing to respond to conventional treatments.

ABSTRACT

IMPACT OF PLATELETS CONCENTRATE AND KERATINOCYTE SUSPENSION ON WOUND HEALING – A PROSPECTIVE RANDOMIZED TRIAL

W. Raffoul, MD¹, S. Guerid, MD¹, D. Hayoz, MD², M. Benathan, PhD³

¹Plastic and Reconstructive Surgery Service, CHUV, Rue du Bugnon 46, 1011 Lausanne, Vaud, Switzerland

²Internal Medicine Department, CHUV, Rue du Bugnon 46, 1011 Lausanne, Vaud, Switzerland

³Skin graft laboratory, Av. de Beaumont 29, 1011 Lausanne, Vaud, Switzerland

BACKGROUND:

Despite the progress achieved in the last decade, wound healing remains a difficult subject to which modern medicine does not always have an efficient answer. Currently, it has been well established that platelets play a fundamental role in wound healing, primarily through the formation of platelet clots and the release of many growth factors and cytokines.

DESIGN:

Prospective randomized trial

METHODS:

- I) The platelet concentrate collected from 10 volunteers using the RegenKit[®] filter system (Regen Lab, Switzerland) was, at first, qualitatively evaluated in vitro.
- II) The study included thin skin graft samples collected from 45 patients between June 2005 and March 2006 at CHUV, Vaud, Switzerland. After randomization, patients were treated either with a standard dressing or using platelet concentrate with or without a suspension of keratinocytes. Healing of the region where the skin graft was first collected was compared across the 3 groups, looking at the duration and quality of the healing process as well as the inflicted pain.

RESULTS:

- I) Blood centrifugation using the RegenKit[®] system (Regen Lab, Switzerland) yielded four times the concentration of platelets compared to standard methods of platelet collection, stability over time and neutrality with regards to coagulation parameters.
- II) Reduction of healing time from 12 days for the control group to 7 days for the APG group and down to 5 days for the APG+K group with significant decrease of inflicted pain.

Fig 1: Duration in days until complete cicatrisation



Fig 2: Epithelialization quality at day 5


Fig 3: Pain evaluation at day 5

CONCLUSIONS:

The results of this study confirm our speculations on the major role of the enriched platelet clot in stimulating wound healing. This also supports the possibility of using this concentrate to carry a suspension of autologous keratinocytes. This latter combination proves to be the most effective way to heal a wound.

This study has also revealed the beneficial effects this product has on pain, inflammatory reactions and the final appearance of the scar. It is interesting to note that the procedure used opens the path towards a more efficient treatment of chronic wounds such as ulcers, eschars, diabetic neuropathies or necrotic angiodermatitis.

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ABSTRACT

THE ROLE OF THE AUTOLOGOUS PLATELET-DERIVED GROWTH FACTOR IN THE MANAGEMENT OF DECUBITUS ULCER.

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ABSTRACT


RegenPRP ENHANCES PROLIFERATION OF HUMAN EPIDERMAL, MESODERMAL AND DERMAL CELLS EX VIVO.

D.F. du Toit¹, W.J. Kleintjes., E.J Mazyala, L du Toit, B.J. Page.

*Division of Anatomy and Histology, Faculty of Health Sciences, University of Stellenbosch.¹
Diabetes Discovery Platform., MRC, South Africa.*

The regenerative outcomes of autologous fibroblast, keratinocyte and skeletal myoblasts and use of GF, have been reported in open ended phase-1 trials. Tissue culture (TC) studies, show that Platelet-Rich Plasma (PRP) can enhance gene expression of matrix molecules such as collagen. PRP is capable of stimulating fibroblast-generated collagen in vitro, due to the trophic platelet growth factor releasate (containing cytokines). Our study was aimed at ex vivo human epidermal, dermal, and mesodermal cell proliferation by RegenPRP, standard and serum free mediums (DMEM). TC were set up by explant technology and static, air-lifting, initially on adhesion molecule covered Petri-dishes, cultured at 37 degrees, 5% CO₂ and 95% air. TC between 7 and 14 days (fibroblasts, keratinocytes, and myoblasts) were replated on RegenPRP or collagen coated plates and exposed to RegenPRP (10,20,40%) and FCS (10,20%) and serum free medium. 100 specimens were equally divided into 415 explants, 3D culturing in fibrin. Results showed that RegenPRP added to DMEM was a far more potent mitogen for the cell line studies, compared to the proliferative capacity of conventional mediums. Also more potent than our historical controls with bFGF. RegenPRP strongly stimulated cell proliferation in the presence of monolayers (keratinocytes, fibroblasts and skeletal myoblasts). Further studies are in progress to study PRP affect on cell line differentiation and phenotype. In conclusion, RegenPRP provides potent trophic mitogens in vitro, devoid in conventional

CARDIO CHIRURGIA

ABSTRACT

THE USE OF AUTOLOGOUS PLATELET GEL TO TREAT DIFFICULT-TO-HEAL WOUNDS: A PILOT STUDY.

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
Tab 1

	treated	control
Group 1: Dehiscent sternal wounds		
Number	10	12
Mean age +/- SD, years	64+/-8	66+/- 5
Male/Female ratio	6/4	8/4
Primary disease, number (%)		
Ischemic disease	8 (80)	9 (75)
Valvular defect	2 (20)	3 (25)
Coexisting disease, number (%)		
Diabetes	5 (50)	7 (58)
Endocarditis	1 (10)	NA
Mediastinitis	1 (10)	NA
None	3 (30)	5 (41)
Group 2: Skin necrotic ulcers		
Number	17	14
Mean age +/- SD, years	61+/- 18	63 +/- 16
Male/Female ratio	8/9	5/9
Primary disease, number (%)		
Venous ulcers	3 (17)	5 (35)
Post-traumatic ulcers	5 (29)	3 (21)
Paraplegic ulcers	3 (17)	2 (14)
Arteriopathic ulcers	1 (6)	2 (14)
Leg- foot pressure ulcers	5 (29)	2 (14)
Coexisting disease, number (%)		
Diabetes	2 (11)	2 (14)
Lymphoma	1 (6)	NA
Haemolytic anemia	1 (6)	NA
Multiplemyeloma	NA	1 (7)
None	10 (59)	9 (64)

CHIRURGIA ODONTOSTOMATOLOGICA

ABSTRACT

BONE REGENERATION OF DENTAL IMPLANT DEHISCENCE DEFECTS USING A CULTURED PERIOSTEUM MEMBRANE

Daiki Mizuno, Hideaki Kagami, Hirokazu Mizuno, Junji Mase, Kazutada Usami, Minoru Ueda, et al.

Hideaki Kagami, DDS, PhD - Division of Stem Cell Engineering - The Institute of Medical Science - The University of Tokyo

OBJECTIVES:

This study aimed to demonstrate the feasibility of a cultured periosteum (CP) membrane for use in guided bone regeneration at sites of implant dehiscence.

Material and methods: Four healthy beagle dogs were used in this study. Implant dehiscence defects (4 × 4 × 3mm) were surgically created at mandibular premolar sites where premolars had been extracted 3 months back. Dental implants (3.75mm in diameter and 7mm in length) with machined surfaces were placed into the defect sites (14 implants in total). Each dehiscence defective implant was randomly assigned to one of the following two groups:

- (1) PRP gel without cells (control) or
- (2) a periosteum membrane cultured on PRP gel (experimental).

Dogs were killed 12 weeks after operation and nondecalcified histological sections were made for histomorphometric analyses including percent linear bone fill (LF) and bone-to-implant contact (BIC).


RESULTS:

Bone regeneration in the treatment group with a CP membrane was significantly greater than that in the control group and was confirmed by LF analysis. LF values in the experimental and the control groups were 72.36 ± 3.14% and 37.03 ± 4.63%, respectively (P<0.05). The BIC values in both groups were not significantly different from each other.

The BIC values in the experimental and the control groups were 40.76 ± 10.30% and 30.58 ± 9.69%, respectively (P=0.25) and were similar to native bone.

Conclusion: This study demonstrated the feasibility of a CP membrane to regenerate bone at implant dehiscence defect.

MEDICINA E CHIRURGIA ESTETICA

ABSTRACT

NEW INJECTION TECHNIQUE BASED ON PATIENT CELL BIOSTIMULATION FOR AUGMENTATION AND REGENERATION WITH REGENACR PROCESS.

Donald du Toit: Research Fellow , Photo-light Biology Laboratory , Faculty of Health Sciences ,University of Stellenbosch and Medical Research Council, Parow, South Africa.

INTRODUCTION:

For predicted outcome with facial rejuvenation techniques, it is imperative to perform prior advanced skin analysis , and understand histologic diagnostic characteristics before regeneration either by the biological or light-based routes. Most often underlying conditions, texture , colour and secretions are neglected by the therapist. This explains why non-ablative thermolysis , fractional laser and resurfacing by IPL or RF are associated with negative factors and outcomes. These include slow gradual repair, variability, need for repeated treatments, scarification, post inflammatory hyperpigmentation, hypopigmentation, non-response to ameliorate horizontal or vertical crepe wrinkles of the face and décolleté. Full-facial rejuvenation at one session, is now possible by the biologic route utilizing biostimulation and autologous platelet-derived growth factors, based on open ended phase-1 studies.

METHODS:

The process is as follows. Skin analysis of the client destined for facial rejuvenation by the biological route, is performed. Venapuncture is performed and 16 ml blood is collected in special anticoagulated tubes supplied by RegenLab, Mollen, Switzerland. Centrifugation follows and PRP and PPP is separated. Prior to rejuvenation the PRP and thrombocytes therein, are activated with calcium chloride. The enriched PRP is used for dermal regeneration by multiple cutaneous and dermal injections, alike mesotherapy, delivering the platelet derived growth factors such as PDGF, TFG, IGF, EGF and VEGF into the papillary dermis and hypodermis. A 30 gauge needle, mounted on a syringe supplied in the RegenKit, is used for autologous cellular regeneration (ACR). The growth factors stimulate resident stem cells and other supporting cells in the ECM. Rejuvenation of the epidermis, dermis with improvement in texture, colour and secretions is the end-point.

RESULTS:

Ex vivo cell biology studies and tissue culture, using Regenlab PRP demonstrate enhanced fibroblast and keratinocyte proliferation. Cell signalling and crawling is stimulated by PRP culture technologies. Clinical studies , and interim results in 10 persons, varying from 40-50 years, of both gender, show safety, with minimal down time of 12-24 hours. Localised swelling subsides, rapidly with the application of ice packs. Improvement of the skin texture, tone and colour is evident within 14 days and is ongoing without need for light-therapy. The response can be enhanced by the adjunct use of light emitting diodes. Subjective patient feedback is encouraging. Skin changes in the glabellar, forehead and nasolabial folds are visible within 14 days by serial photography. The longest follow –up is 18 months. A small percentage may need directed re-injection for deep lines.

CONCLUSION:

Autologous PRP, as utilised in the RegenlabACR biological process, contains potent trophic factors relevant to facial regeneration, demonstrable ex vivo in tissue culture and in the clinic, portraying minimal side effects.

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→ DR. VÍCTOR GARCÍA GIMÉNEZ

E-mail:
director@cmadieuropa.com

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- Sanchez M, Azofra J, Arntua E, Andia I, Padilla S, Bantisteben J, Mujika I. Plasma Rich In Growth Factors To Treat An Articular Cartilage Avulsion: A Case Report. Arthroscopic Surgery Unit, USP-La Esperanza Clinic, Vitoria-Gasteiz, Basque Country, SPAIN. *Med Sci Sports Exerc* 2003;35(10):1648-1652
- Alex Dugrillon and Harald Küller. Topical Application Of Platelets For Improved Wound Healing. Institute of Translational Medicine and Immunology, Faculty of Clinical Medicine Mannheim, University of Heidelberg, German Red Cross Blood Services of Baden-Württemberg, Mannheim, Germany. *Blood Therapies in Medicine* 2002;3(1)
- Arntua E. The Use Of Plasma-Rich In Growth Factors (Prpf) In Oral Surgery. *Pract Proced Aesthet Dent*. 2001 Aug;13(8):487-93; quiz: 487-93/*Pract Proced Aesthet Dent* 2001;13:487-93
- Arntua E. Expansión De Cresta Con Osteotomías: Estado Actual. Utilización Del Plasma Rico En Factores De Crecimiento (PRGF). *Rev Esp Chirg Oral y Maxilofacial* 2001;23:1-5
- Arntua E. Factores De Crecimiento Plasmáticos. Una Revolución Terapéutica. *Ideas y Trabajos Odontostomatológicos* 2001;2(2):90-94
- Welsh, W.J. Platelet Gel Usage In Cosmetic Surgery. *American Academy of Cosmetic Surgery*, 2000;27
- Welsh, W.J. Autologous Platelet Gel. Clinical Function and Usage In Plastic Surgery. *Cosmetic Dermatology* 2000;13-18
- Croce MA, D'Yna K, Boraldi F, Guaglio D Jr, Cetta G, Tiozzo R and Pasquali Ronchetti I. Hyaluronan affects protein and collagen synthesis by in vitro human skin fibroblasts. *Tissue Cell* 2001;33(4):326-31
- Wang Q, Lu K, Yang L. Effects of hyaluronic acid-stimulating factor on viability and collagen synthesis of fibroblasts
- Zhonghua Zheng Xing Shao Shang Wai Ke Za Zhi. 1999 Mar;15(2):99-101
- Isenard N, Legault JM, Renard G, Robert L. Effect of hyaluronan on MMP expression and activation. *Cell BioInt* 2001;25(8):736-9
- Lesley J, Hascall VC, Tammi M, Hyman R. Hyaluronan binding by cell surface CD44. *J Biol Chem* 2000;1;275(35):25067-75
- Raffitt DM, Ogston N, Jugdeochsingh R, Cheung HF, Evans BA, Thompson RR, Powell JJ, Hampson GN. Orthosilicic acid stimulates collagen type 1 synthesis and osteoblastic differentiation in human osteoblast-like cells in vitro. *Bone* 2003;32(2):127-35
- Cattemme MR, Vanden Berghhe DA. Supplementation of calves with stabilized orthosilicic acid. Effect on the B1, Ca, Mg, and P concentrations in serum and the collagen concentration in skin and cartilage. *Biol Trace Elem Res* 1997;56(2):153-65
- Carlisle EM. In vivo requirement for silicon in articular cartilage and connective tissue formation in the chick. *J Nutr* 1976;106(4):478-84
- Shi X, Ding M, Chen F, Wang L, Rojanasakul Y, Vallyathan V, Castaneda V. Reactive oxygen species and molecular mechanism of silica-induced lung injury. *J Environ Pathol Toxicol Oncol* 2001;20(Suppl 1):85-93
- Baroni T, Bodo M, D'Alessandro A, Conte C, Calvitti M, Muzi G, Lumina A, Bellocchio S, Abbruti G. Silica and its antagonistic effects on transforming growth factor-beta in lung fibroblast extracellular matrix production. *J Invest Med* 2001;49(2):146-56
- Hubbard AK, Thibodeau M, Giardino C. Cellular and molecular mechanisms regulating silica-induced adhesion molecule expression in mice. *J Environ Pathol Toxicol Oncol* 2001;20(Suppl 1):45-51
- Grossman R. The role of dimethylaminoethanol in cosmetic dermatology. *Am J Clin Dermatol* 2006;6(1):39-47
- Kanchwala BK, Holloway L, Bucky LP. Reliable soft tissue augmentation: a clinical comparison of injectable soft-tissue fillers for facial-volume augmentation. *Ann Plast Surg* 2006;55(1):30-5; discussion 35
- Mamur EB, Phelps R, Goldberg DJ. Clinical, histologic and electron microscopic findings after injection of a calcium hydroxylapatite filler. *J Cosmet Laser Ther* 2004;6(4):223-6
- Flaherty P. Radience. *Facial Plast Surg*. 2004;20(2):165-9
- Comite SL, Liu JF, Balasubramanian S, Christian MA. Treatment of HIV-associated facial lipotrophy with Radience FN (Radless). *Dermatol Online J* 2004;10(2):2
- Tzikas TL. Evaluation of the Radience FN soft tissue filler for facial soft tissue augmentation. *Arch Facial Plast Surg* 2004;6(4):234-9
- Lemperle G, Morham V, Charrier U. Human histology and persistence of various injectable filler substances for soft tissue augmentation. *Aesthetic Plast Surg* 2003;27(3):354-66; discussion 367
- Chodorowska G, Tomczyk M, Glowacka A. Basic-fibroblast growth factor (b-FGF)—its biological role in physiologic and pathologic conditions. *Ann Univ Mariae Curie Skłodowska [Med]*. 2004;59(1):286-91

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- Fournier PF. *Lipoesculture: the syringe technique*; 1991
- Peck H, Peck S. A concept of facial aesthetics. *Anglen Orthot* 1970;40:284-317
- Deweese DD, Saunders WH, Shuller DE, et al. *Otolaryngology, head and neck surgery*. St. Luis, Mosby, 1988
- Fournier PF. Body sculpturing through syringe liposuction and autologous fat re-injection. *Samuel Rolf International* 1987
- Daly CH, Odland GF. Age-related changes in the mechanical properties of human skin. *J Invest Dermatol* 1979; 73(1):84-87
- Toriumi DM, East CA, Larrabee WF. Osteoinductive biomaterials for medical implantation. *J. Long-term Effects Med Implant* 1991;1(1):53
- Agris J. Autologous fat transplantation and suction assisted lipectomy. *A Clinical Atlas*, 2a ed. Houston: Eclectic Enterprises, 1987
- Fournier PF. *Fattransfer*. En: Parish LC, Lask GP (eds). *Aesthetic Dermatology*. Nueva York: McGraw Hill, 1992;2:275
- Hernández I, Rozan G, Dávila M. Reingeniería de tejidos. Plasma rico en plaquetas como inductor de reparación en la panculopatía edematofibroesclerótica. *American Journal of Medicine and Cosmetic Surgery* 2005;7(2):18-24
- Bircoll M. A Nine-year experience with autologous fat transplantation. *Am J Cosmet Surg* 1992;9:55-59

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