



Identification of *Anisakidae* Larvae at the species level by multiplex PCR

INDEX

1. Aim and field of application	2
2. Principle of the method	2
3. References	2
4. Definitions	3
5. Devices/instruments	3
6. Reagents and chemicals	4
7. Procedure	
7.1 Sample preparation	5
7.2 Method	5
8. Results	8
9. Characteristics of the method	8
10. Safety measures	8

1. Aim and field of application

To identify the species or genotype of single *Anisakidae* larvae preserved in ethanol by a multiplex PCR analysis. This method can be applied to larvae, or their portions, collected from human biopsies or from tissues of animal origin.

2. Principle of the method

The PCR is a molecular biology technique that allows for the amplification of specific nucleic acid fragments, which initial and terminal nucleotide sequences are known (oligonucleotide pair). If a species (or genotype) has its own characteristic DNA portion, due to its composition and/or dimension, it is possible to choose an oligonucleotide pair allowing for its amplification. The PCR amplification is characterized by a high sensitivity and specificity. Multiplex-PCR is an application of standard PCR where two or more pairs of primers are used at the same time. In such a way, it is possible to amplify multiple sequences at the same time during the same reaction.

The larval stages of nematodes of the *Anisakidae* family parasitize fishes, cephalopods and shrimps and the *Anisakidae* adult worms parasitize fishes, marine birds and mammals. The larvae of the genus *Anisakis* and *Pseudoterranova* are responsible, if ingested by humans, of the disease known as anisakiasis. Larvae, obtained from human patients or infected hosts, can be morphologically identified, however larvae of the genus *Anisakis* and *Pseudoterranova* cannot be distinguished.

Multiplex PCR molecular method, based on the different size of the PCR fragments amplified from a portion of the ITS (Internal Transcribed Spacer) region, allowed to distinguish single larvae of *Anisakis* spp. from that of *Pseudoterranova*, *Contracaecum* and *Hysterothylacium* spp., and, within *Anisakis* spp., between *A. pegreffii*, *A. simplex* s.l. (including the hybrid genotype *A. simplex/A. pegreffii*), *A. physeteris* (including *A. brevispiculata* and *A. paggiae*) and *A. typica* at the species level.

The sizes of the fragments produced by the amplification with oligonucleotide pairs specific for each member of the *Anisakidae* family, are shown in Table A.

Table A – Expected size of the ITS fragments (base pairs) for each species obtained with the reported primers.

Primer 5'	Primer 3'	<i>A. pegreffii</i>	<i>A. simplex</i> s.l. e <i>A. simplex/pegreffii</i> ibrido	<i>A. physeteris</i> ^a	<i>A. typica</i>	<i>Contracaecum</i> <i>oscultatum</i>	<i>C. rudolphii</i> (A, B, C)	<i>Pseudoterranova</i> spp.	<i>Hysterothylacium</i> <i>aduncum</i>
ASPf	RevB	588bp	588bp						
APE1	"	672bp							
APyf	"			143bp					
ATf	"				427bp				
COf	"					799bp			
CRf	"						307bp		
PDf	"							370bp	
HAf	"								991bp

^a can amplify also *Anisakis brevispiculata* and *A. paggiae*

3. References

La Rosa G., D'Amelio S e Pozio E. (2006) "Molecular Identification of nematode Worms From Seafood (*Anisakis* spp. and *Pseudoterranova* spp.) and Meat (*Trichinella* spp.), in *Food-Borne Pathogens. Methods and Protocols*. Edited by (Adley CC., ed.) Human Press, Totowa, NJ, pp. 217-234

Murrel KD. (2002) Fishborne zoonotic parasites: epidemiology, detection and elimination, in *Safety and Quality Issues in Fish Processing* (Bremmer, HA., ed.) CRC Press, Boca Raton, FL, pp.114-141.

Zhu XQ, D'Amelio S, Palm HW, Paggi L, George-Nascimento M, Gasser RB. (2002) SSCP-based identification of members within the *Pseudoterranova decipiens* complex (Nematoda: Ascaridoidea: Anisakidae) using genetic markers in the internal transcribed spacers of ribosomal DNA. *Parasitology*. 124, pp. 615-23.

D'Amelio S, Mathiopoulou KD, Santos CP, Pugachev ON, Webb SC, Picanço M, Paggi L. (2000) Genetic markers in ribosomal DNA for the identification of members of the genus *Anisakis* (Nematoda: ascaridoidea) defined by polymerase-chain-reaction-based restriction fragment length polymorphism. *Int J Parasitol.* 30, pp. 223-226.

Zhu X, Gasser RB, Podolska M, Chilton NB. (1998) Characterisation of anisakid nematodes with zoonotic potential by nuclear ribosomal DNA sequences. *Int J Parasitol.* 28, pp. 1911-21.

Kijewska A, Rokicki J, Sitko J, Wegrzyn G. (2002) Ascaridoidea: a simple DNA assay for identification of 11 species infecting marine and freshwater fish, mammals, and fish-eating birds. *Exp Parasitol.* 101, pp.35-39.

Mattiucci S, Nascetti G. (2006) Molecular systematics, phylogeny and ecology of anisakid nematodes of the genus *Anisakis* Dujardin, 1845: an update. *Parasite.* 13, pp. 99-113.

Farjallah S, Busi M, Mahjoub MO, Slimane BB, Paggi L, Said K, D'Amelio S. (2008) Molecular characterization of larval anisakid nematodes from marine fishes off the Moroccan and Mauritanian coasts. *Parasitol Int.* 57, pp.430-6.

ISO 22174: 2005. Microbiology of food and animal feeding stuffs – Polymerase chain reaction (PCR) for the detection of food-borne pathogens.

ISO/FDI 20837:2006(E). Microbiology of food and animal feeding stuffs – Polymerase chain reaction (PCR) for the detection of food-borne pathogens - Requirements for sample preparation for qualitative detection

ISO/FDI 20838:2006(E). Microbiology of food and animal feeding stuffs – Polymerase chain reaction (PCR) for the detection of food-borne pathogens - Requirements for amplification and detection for qualitative methods

Promega Technical Bulletin, Part# TB307 Rev 5/06. Tissue and hair extraction kit (for use with DNA IQ™) Protocol. Instruction for use of product DC6740

Promega Technical Bulletin, Part# TB297 Rev 4/06. DNA IQ™ system – Database Protocol. Instruction for use of products DC6700 and DC6701

Reinecke, S.N., Morgan, R.D. (1991) BfaI, a new *MaeI* isoschizomer from *Bacteroides fragilis*, recognizes the sequence 5' C decreases TAG 3'. *Nucleic Acids Res.* 19, pp. 1152

Ono, A., Matsuo, Y., Matsuda, A., Ueda, T. (1993) Nucleosides and nucleotides. CXIX. Inhibition of DNA-cytosine methylase HhaI by a self-complementary oligonucleotide containing 5-fluorocytosine. *Biol. Pharm. Bull.* 16, pp. 529-533

Sato, S., Hutchison, C.A. III, Harris, J.I. (1977) A thermostable sequence-specific endonuclease from *Thermus aquaticus*. *Proc. Natl. Acad. Sci. USA* 74, pp. 542-546

Abollo E, Paggi L, Pascual S, D'Amelio S. (2003) Occurrence of recombinant genotypes of *Anisakis simplex* s.s. and *Anisakis pegreffii* (Nematoda: Anisakidae) in an area of sympatry. *Infect Genet Evol.* 3, pp. 175-181.

Umehara A, Kawakami Y, Araki J, Uchida A. (2008) Multiplex PCR for the identification of *Anisakis simplex* sensu stricto, *Anisakis pegreffii* and the other anisakid nematodes. *Parasitol Int.* 57, pp. 49-53.

4. Definitions

ITS (Internal Transcribed Spacer 1), interspaced sequence of the nuclear ribosomal gene including the ITS-1, the 5.8S gene, ITS-2 gene with the addition of 70 bp of the 28S gene

Oligonucleotide, short sequence (15/30 nucleotide bases) used to amplify a DNA specific fragment

Multiplex Primers Mix, mix of 9 oligonucleotide base pairs amplifying ITS from each species

Reference larvae, larvae of *Anisakis pegreffii* preserved in ethanol

Reference DNA, DNA extracted from reference larvae, supplied

Positive control for the DNA extraction, a reference larva analysed in the same working session of test samples, to verify the efficacy of the DNA extraction session

DNA/larva, DNA extracted from a single larva

Positive control for the amplification, a reference DNA; this control is used in the amplification session to verify the efficacy of the PCR

Negative control for the amplification, reagent grade water; this control is used in the amplification session to verify the absence of contamination in the PCR reaction

PCR, Polymerase Chain Reaction.

The definitions and terminology used in the UNI EN ISO 22174 standard are applied in the present method.

5. Devices/instruments

- 5.1 Stereo microscope, magnification 60÷100x
- 5.2 Bench top refrigerated centrifuge for 1.5 mL tubes, minimum 10,000xg
- 5.3 Freezer ≤-15°C
- 5.4 Thermoblock with vibration, temperature range 25÷100°C
- 5.5 Magnetic separation stand
- 5.6 PCR thermocycler
- 5.7 Refrigerator, temperature range 1÷8°C
- 5.8 Qiaxcel, capillary electrophoresis system
- 5.9 Horizontal electrophoretic apparatus
- 5.10 Analytical balance, readability 0.1 g
- 5.11 UV transilluminator
- 5.12 Digital imaging system
- 5.13 Adjustable volume pipettes, volume range: 1-10 µL, 2-20 µL, 20-100 µL, 50-200 µL, 200-1000 µL
- 5.14 Reagent grade water system production
- 5.15 Vortex
- 5.16 Orbital shaker

6. Reagents and chemicals

- 6.1 **Incubation buffer.** Commercial solution: DNA IQ™ System kit, Promega, code DC6701 or DC6700. Once prepared, label the solution with "IB+". Store PK and DTT and their aliquots at -20°C, whereas store the incubation buffer at +4°C.
- 6.2 **Lysis buffer.** Commercial solution: Tissue and Hair Extraction Kit, Promega, code DC6740. Once prepared, label the solution with "LB+". Store DTT and its aliquots at -20°C, whereas store the lysis buffer at room temperature.
- 6.3 **Paramagnetic resin.** Commercial suspension: DNA IQ™ System kit, Promega, code DC6701 or DC6700. Store at room temperature.
- 6.4 **Washing buffer.** Commercial solution: Tissue and Hair Extraction Kit, Promega, code DC6740. Once prepared, label the solution with "WB+". Store at room temperature
- 6.5 **Eluting buffer.** Commercial solution: Tissue and Hair Extraction Kit, Promega, code DC6740. Store at room temperature
- 6.6 **PCR master mix Hotstart.** Commercial solution. Store according to the manufacturer's recommendations.
- 6.7 **Oligonucleotides.** Commercial preparation (*Table B*); the lyophilized products is reconstituted with TE 0.1x, according to the manufacturer's recommendations, at a concentration of 100 pmol/µL; the

lyophilized product can be stored frozen for up to 20 years; the reconstituted product can be stored frozen up to 10years.

- 6.8 10X multiplex Primers Mix.** The oligonucleotides mixture (6.7) used for the multiplex PCR; the mixture is obtained combining a defined volume of the oligonucleotides reported in *Table B* in TE 0.1x (6.16); the final concentration corresponds to 3 pmol/μL, except for APEf, at 4 pmoli/μL, and RevB, at 6 pmoli/μL. 100μL aliquots are prepared and stored frozen up to 10years.

Table B – Oligonucleotides present in the 10X Multiplex Primers' Mix (6.7 e 6.8), their codes and amplified nucleotide sequences.

Oligonucleotide sequence	Code	Target sequence
5'- TTGCAATCACTTCTCTCAGATTG -3'	ATf	ITS
5'- GGCTGGTTGATGAACTGTTG -3'	APyf	ITS
5'- GACATTGTTATTTTCATTGTATGTGTTGAAAATG -3'	ASPf	ITS
5'- GAGCAGCAGCTTAAGGCAGAGGC -3'	APE1	ITS
5'- CGAGTACTTTTTATGGTCGTGAAGT -3'	PDf	ITS
5'- TGATATGCTTGAAAGGCAGG -3'	COf	ITS
5'- CGACAAGCAGTGTCCCTTTG -3'	CRf	ITS
5'- GCCTTCCATATGCGCGTATA -3'	HAf	ITS
5'- GCCGGATCCGAATCCTGGTTAGTTTCTTTTCC-3'	RevB T	ITS

- 6.9 Loading buffer.** Commercial product allowing DNA molecule electrophoresis to be performed. Store according to the manufacturer's recommendations.
- 6.10 Agarose.** Commercial product suitable for performing DNA molecule electrophoresis. Store at room temperature according to the manufacturer's recommendations.
- 6.11 TAE solution 50x.** Commercial product (2M Tris-acetate, 50mM EDTA, pH 8.2–8.4 at 25°C). Store according to the manufacturer's recommendations.
- 6.12 TAE solution 1x.** 1000 mL preparation: take 20 mL of the 50x solution and bring to 1000 mL with water. Store at room temperature for up to 1 month.
- 6.13 DNA intercalating agent.** Commercial product able to fit in the DNA double helix, used to view the amplification product on agarose gel. Store according to the manufacturer's recommendations..
- 6.14 L50.** Commercial product containing markers for DNA molecular weight multiple of 50 bp. All commercial products containing molecules of multiples of 50 bp within the 50-1000 bp range can be used. Store refrigerated according to manufacturer's recommendations.
- 6.15 TE 1x solution.** Commercial product 10mM Tris-HCl (pH 8,0), 1mM EDTA-Na₂, pH 7.9–8.1 at 25°C. Store refrigerated according to the manufacturer's recommendations.
- 6.16 TE 0.1x solution.** TE 1x diluted solution. For 100 mL preparation: take 10 mL from 1x solution and add 90 mL water. Filter with 0.22 μm filters and prepare aliquots of 10 mL. Store frozen for up to 24 months.
- 6.17 QIAxcel DNA High Resolution kit.** Commercial product to be used with QIAxcel (5.8). It contains high-resolution gel cartridge and buffers for sample preparation and analysis. Store according to manufacturer's recommendations.
- 6.18 Alignment marker 15-600 bp** (Alignment marker with 15 bp and 600 bp fragments).Commercial product to be used with QIAxcel (5.8). Store according to manufacturer's recommendations.
- 6.19 DNA size marker.** Commercial product to be used with QIAxcel (5.8). Store according to manufacturer's recommendations.

6.20 Reagent grade water or Milli-Q.

6.21 Reference Larvae. Larvae or larvae portion of *A. pegreffi*, stored in ethanol (95-99%). Store frozen for up to 10 years.

6.22 Reference DNA. Genomic DNA purified from a reference larvae. Store frozen for up to 10 years.

7. PROCEDURE**7.1 Sample preparation**

- Test samples are inspected to verify the presence of larvae and the preservation conditions. If the condition are not suitable, the test is not performed.
- The ethanol containing the larvae is transferred into a Petri dish and observed under the stereo microscope. A maximum of three larvae are collected and placed in 1.5 mL conical tubes, one larva in each tube. Excess ethanol is removed and the minimum volume is left.
- Spin tubes containing larvae at maximum speed for a few seconds.
- Store the tubes frozen. Under these conditions, larvae can be stored for the DNA extraction for up to 10 years.

7.2 Method**7.2.1 DNA extraction from single larva**

- If not otherwise specified, the procedure is carried out at room temperature.
 - Each working session requires the DNA extraction of a reference larva identified as “positive control for the extraction”.
 - Before starting the procedure, prepare a sufficient volume of the IB+ (6.1) and LB+ (6.2) solutions according to the manufacturer’s recommendations.
- a) Centrifuge the tubes containing the larvae to be identified at maximum speed for a few seconds.
 - b) Add 100 µL of IB+ (6.1).
 - c) Incubate at 55°C for 30-60 min in the thermoblock. During incubation, shake at 1,400 vibrations/min.
 - d) Centrifuge, as in point “a”.
 - e) Add 200 µL of LB+ (6.2).
 - f) Add 10 µL of paramagnetic resin (6.3) after resuspending it by vortexing.
 - g) Incubate for 5-10 min at 25°C in the thermoblock. During incubation, shake at 1,400 vibrations/min.
 - h) Place the tubes in the magnetic separation stand and wait for 30-60 sec, so that the paramagnetic resin particles are blocked by the magnetic field on the tube wall.
 - i) Discard the liquid phase by aspirating, avoiding the dislodging of the resin particles.
 - j) Add 100 µL of LB+ (6.2) and resuspend the resin particles by vortexing.
 - k) Place the tubes in the magnetic separation stand, as in point “h”.
 - l) Discard the liquid phase by aspirating.
 - m) Add 100 µL of WB+ 1x (6.4) and resuspend the resin particles by vortexing.
 - n) Place the tubes in the magnetic separation stand, as in point “h”.
 - o) Discard the liquid phase by aspirating.
 - p) Repeat the washing step, from “m” to “o”, with WB+ (6.4) three times.
 - q) After the last washing, leave the tubes open to let the resin particles dry for 15-20 minutes.
 - r) Add 100 µL of the eluting buffer (6.5) and gently resuspend the resin particles, do not vortex.
 - s) Incubate at 65°C for 5 minutes. During incubation, shake at 1,400 vibrations/min.
 - t) Place the tubes in the magnetic separation stand, as in point “h”.

- u) Collect the liquid phase (about 90-100 µL) and transfer it to a 1.5 mL tube.
- v) The resulting extract is defined as "DNA/larva" and stored frozen. Under these conditions, it can be stored for up to 10 years.

7.2.2 Multiplex PCR amplification

- Unless otherwise clearly stated, store tubes on ice; use tips with barrier and wear disposable gloves.
 - At each working session, use a positive and a negative amplification control. Use reference DNA (6.22) as positive control and water (6.20) as negative control.
- a) Thaw DNA/larva, 2x PCR MasterMix (6.6), 10X multiplex Primers Mix (6.8), and positive amplification controls (6.22).
 - b) Mark with a progressive number an adequate number of 0.2 mL PCR tubes.
 - c) Prepare an adequate cumulative volume of the amplification mix. Evaluate the volume on the basis of a single sample amplification mix (*Table C*) and of the total number of samples added of three reactions (one for the positive amplification control, one for the negative one and an extra reaction).

Table C – single sample amplification mix: components and volumes

2x PCR MasterMix (6.6)	25 µL
10x Multiplex PCR Primers Mix (6.8)	5 µL
H ₂ O (6.20)	22 µL
Total	48 µL

- d) Mix the amplification mix by vortexing and centrifuge at maximum speed for a few sec.
- e) Transfer 48 µL of the cumulative amplification mix to each PCR tube (point "b").
- f) Add 2 µL of the larval DNA to be tested to each tube.
- g) Close the tubes, mix by vortexing and centrifuge at maximum speed for a few sec.
- h) Start the amplifying cycle (*Table D*) on the thermocycler device; wait for the temperature to reach 95°C and insert the tubes in the thermoblock by pausing the instrument.

Table D – amplification cycles

Pre-denaturation #	2 min/95°C
Amplification	30 s/95°C 30 s/52°C 45 s/72°C
Number of cycles	30
Final extension	7 min/72°C

pre-denaturation length may vary according to the Master Mix manufacturer's recommendations

- i) At the end of the amplification phase, centrifuge the tubes at maximum speed for a few sec.
- l) Add loading buffer (6.9) according to the manufacturer instructions.
- m) Vortex and centrifuge the tubes at maximum speed for a few sec.
- n) Keep tubes on ice or refrigerated until starting electrophoresis.

7.2.3 Result display by gel electrophoresis

- a) Switch on the Qiaxcel instrument (5.8) and the relative Qiaxcel ScreenGel management software on the PC;
- b) access the "Process Profile" panel; indicate under "Cartridge Type" the option "DNA HighRes"; indicate the desired profile and Experiment Directory;
- c) move the tube tray to the "Access Position" by selecting the "Load Position" item from the "Status Information" panel;

- d) insert in MARKER1 position the 12 tubes containing at least 10 µL of the chosen "Alignment Marker" (6.18); then return the tray to the initial position by selecting the "Park Position" from the "Status Information" panel;
- e) starting from the "A" row, position the samples to analyze (minimum volume 10 µL) in rows of 12 tubes. If the samples to be analyzed are not enough to complete the row of 12, add the required number of tubes containing the QX DNA dilution buffer (minimum volume 10 µL) supplied with the QIAxcel DNA High Resolution kit (6.17)
- f) for each analysis session (which may include up to a maximum of 8 runs of 12 samples), include a tube containing the DNA size marker (6.19)
- g) in "Run Parameters" set the option 0M500 under "Method"; select the runs on the virtual plate in the "Sample Row Selection" side panel;
- h) in "Sample Selection" set the run parameters as follows: "Plate ID": PCR + data "Alignment Marker": 15bp-600bp (6.18). In "Sample Information" enter the names of the samples in the corresponding boxes.
- i) in "Run Check" check that all the selected rows are occupied by tubes containing PCR products to be analyzed or QX DNA dilution buffer (6.17) and that the alignment Marker has been loaded, then select the appropriate boxes; finally select "Run";
- j) visualize the results by selecting the "Absolute migration time" mode from the "Image options" menu and process the data with the "Start analysis" command
- k) scroll through the electropherogram of each sample to check peaks above the highest band of the alignment marker (6.18)
- l) print the results to archive
- m) at the end of the run, close the program and turn off the instrument.

If the instrument 5.8 is out of service for an extended period proceed with agarose gel following the protocol below:

- a) Assemble the electrophoresis apparatus according to the manufacturer's recommendations. For the gel preparation, use a comb suited for the number of samples.
- b) Add 1.6 gr agarose (6.10) in 100 mL TAE 1x (6.12) in a glass beaker.
- c) Gently resuspend the powder by rotation.
- d) Boil the agarose suspension for 30 sec. If the solution is not homogeneous, continue to boil for another 30 sec.
- e) Restore with water the volume lost by boiling.
- f) Allow the agarose solution to cool.
- g) Before it solidifies (at about 47°C), add DNA intercalating agent (6.13) according to the manufacturer's instructions.
- h) Shake gently to dissolve uniformly the DNA intercalating agent and pour the agarose in the gel tray previously prepared (point "a").
- i) Wait for the gel to solidify, which requires at least 30 min.
- j) Place the tray with the gel in the electrophoresis apparatus.
- k) Cover the gel with TAE 1x buffer (6.12) and gently pull out the comb.
- l) The first or last well are loaded with 15 µL of the L50 solution (6.14).
- m) Load in each well 10 µL of the amplification product (point 7.2.2 "i"), respecting the progressive numbering of the tubes (point 7.2.2 "b").
- n) Connect the electrophoresis apparatus with the power supply and set 10 V/cm of gel.
- o) Run the gel for about 30 min or until the fastest dye, contained in the loading buffer (6.9), reaches a distance of 1 cm from the gel border.
- p) After 30 min, switch off the power supply, place the gel under UV illumination and check the band separation. The electrophoresis run is adequate if it is possible to easily distinguish all bands of molecular

weight marker ranging from 250 to 1000 bp. If the separation is incomplete, continue the run.

q) At the end of the run, transfer the gel to the imaging system and print the result.

7.2.4 Result Interpretation by gel electrophoresis

The size of the amplification bands revealed by the electrophoresis is evaluated by comparison with the reference molecular weight L50 (6.14) and with the positive controls of extraction and amplification. The visual evaluation is considered sufficient and adequate (see *Table A*).

The amplification test is considered valid if:

- a) the amplification of the positive control shows an amplification product as in *Table A*;
- b) the amplification of the negative control does not show any amplification product or, eventually, only bands related to unused oligonucleotides and/or primer dimer
- c) the positive control of the extraction product shows an amplification product as in *Table A*;

The species identification is made comparing the size of the amplification fragment(s) produced by the sample(s) with those shown in *Table A*.

In case the sample shows an unexpected band, the identification of the sample at the species level is not possible.

8. Results

The results are expressed as follows:

If the amplification band is comparable with 588bp, the sample is identified as *A. simplex* sl.

If the amplification bands are comparable with the doublet 588bp-672bp, the sample is identified as *A. pegreffi*.

If the amplification band is comparable with 143bp, the sample is identified as *A. physeteris*.

If the amplification band is comparable with 427bp, the sample is identified as *A. typica*.

If the amplification band is comparable with 370bp, the sample is identified as *Pseudoterranova* spp.

If the amplification band is comparable with 799bp, the sample is identified as *Contracaecum osculatum*

If the amplification band is comparable with 307bp, the sample is identified as *Contracaecum rudolphii* (A, B, C).

If the amplification band is comparable with 991bp, the sample is identified as *Hysterotilacium aduncum*.

In case the digestion test was valid but the sample displays a profile of bands not comparable with those reported in *Table A*, the identification at the species level is considered "impossible".

9. Characteristics of the method

This method has been characterised in terms of repeatability and specificity. The results of the validation protocol are used to confirm that the method is suitable for the expected aim and are reported in the validation report, which can be sent upon request.

10. Safety measures

This method has to be carried out only by authorized personnel. The operator should wear individual protection devices during the test performance. For general safety measures, refer to the CDC guidelines.