

**Identification and characterization of
Shiga Toxin (Verocytotoxin)-producing *Escherichia coli* (STEC/VTEC)
by Real Time PCR amplification of the main virulence genes
and the genes associated with the serogroups mainly associated
with severe human infections**

1. Aim and field of application

STEC are *Escherichia coli* strains harboring lysogenic bacteriophages carrying genes encoding the production of Shiga Toxins (Stx) also termed as Verocytotoxins (VTs). The present method is applied to detect by Real Time PCR the presence of the main virulence genes in *E. coli* strain cultures, for their identification as STEC, and of the genes associated with five O serogroup causing the majority of HUS cases in Europe.

The *E. coli* genes comprised in the field of application of this method include:

- 1) stx genes (stx1 and stx2 groups) encoding the Verocytotoxins (Shiga toxins), the main virulence factors of STEC;
- 2) The *eae* gene, encoding a 90KDa protein, the intimin, which is the key factor for the induction of the “attaching and effacing” lesion on the enterocyte, a typical feature of STEC strains belonging to the top-5 serogroups;
- 3) Serogroup-associated genes. The method is aimed at the identification of the STEC serogroups causing the majority of HUS cases in Europe: O157, O26, O111, O103, O145. This is achieved by Real-Time PCR amplification of fragments of genes comprised in the operons encoding the different lipopolysaccharides (LPS) constituting the O antigens or anyhow associated to each serogroup in a unique manner (see the Table below).

This method is applied to DNA purified from isolated bacterial strains as the matrix, for the confirmation of the presence of the virulence genes and serogroup identification.

The target genes and the primers and probes used in this method are the same included in the international standard for the detection of STEC in food ISO/TS 13136: 2012.

2. Procedure

2.1. DNA extraction and purification

One ml of an overnight culture in a non-selective medium suitable for *E. coli* is used for DNA extraction. This step is accomplished by any method in use in the laboratory for the extraction and purification of DNA from food enrichment cultures.

When dealing with cultures from solid media, an isolated colony is suspended in 1 ml of distilled water and used for DNA extraction.

2.2 Real-Time PCR amplification

The protocol is based on the 5' nuclease PCR assay. Considering that Real-Time PCR may use different instruments and probes labelling chemistry, the amplification conditions to be applied may vary depending on the system used. Refer to the instructions supplied with the instrument and kit of choice.

The primers and probes to be used are listed in the tables below. The chemistry of the reporter and quencher fluorophores is not indicated being largely dependent on the Real-Time PCR systems available in each laboratory. The bibliographic references for the primers and probes sequences are indicated in the table captions.

Table: Primer sequences and amplicon sizes.

Target	Forward primer, reverse primer and probe sequences (5'-3')	Amplicon size (bp)	Location within sequence	GenBank accession number
stx1 (Perelle S. <i>et al.</i> Mol Cell Probes 2004 18:185–192)	Fwd- TTTGTACTGTSACAGCWGAAGCYTTA CG Rev- CCCCAGTTCARWGTRAGRTCMACRTC Probe- CTGGATGATCTCAGTGGGCGTTCTTAT GTAA	131	878–906 983–1008 941–971	M16625
stx2 (Perelle S. <i>et al.</i> Mol Cell Probes 2004 18:185–192)	Fwd- TTTGTACTGTSACAGCWGAAGCYTTA CG Rev- CCCCAGTTCARWGTRAGRTCMACRTC Probe- TCGTCAGGCACTGTCTGAAACTGCTC C	128	785–813 785–813 838–864	X07865
eae (Møller Nielsen E. and Thorup Andersen M. J clin Microbiol 2003)	Fwd- CAT TGA TCA GGA TTT TTC TGG TGA TA Rev- CTC ATG CGG AAA TAG CCG TTA Probe- ATAGTCTCGCCAGTATTCGCCACCAAT ACC	102	899-924 1000-979 966-936	Z11541
rfbEO157 (Perelle S. <i>et al.</i> Mol Cell Probes 2004 18:185–192)	Fwd- TTTCACACTTATTGGATGGTCTCAA Rev- CGATGAGTTTATCTGCAAGGTGAT Probe- AGGACCGCAGAGGAAAGAGAGGAATT AAGG	88	348–372 412–435 381–410	AF163329
wbdI O111 (Perelle S. <i>et al.</i> Mol Cell Probes 2004 18:185–192)	Fwd- CGAGGCAACACATTATATAGTGCTTT Rev- TTTTTGAATAGTTATGAACATCTTGTTT AGC Probe- TTGAATCTCCAGATGATCAACATCGT GAA	146	3464–3489 3579–3609 3519–3548	AF078736
§wzxO26 (Perelle S. <i>et al.</i> Mol Cell Probes 2004 18:185–192)	Fwd- CGCGACGGCAGAGAAAATT Rev- AGCAGGCTTTTATATTCTCCAACCTTT Probe- CCCCGTTAAATCAATACTATTTACGA GTTGA	135	5648–5666 5757–5782 5692–5724	AF529080
wzxO145 (USDA, 2012)	Fwd- AAA CTG GGA TTG GAC GTG G Rev- CCC AAA ACT TCT AGG CCC G Probe- TGCTAATTGCAGCCCTTGCACTACGAG GC	132	2600222–2600240 2600106–2600124 2600190–2600162	AP019708
wzxO103 (Perelle S. <i>et al.</i> J Appl Microbiol 2005 98:1162–1168)	Fwd- CAAGGTGATTACGAAAATGCATGT Rev- GAAAAAAGCACCCCGTACTTAT Probe- CATAGCCTGTTGTTTTAT	99	4299–4323 4397–4375 4356–4373	AY532664

2.3. Controls

A positive and two negative controls are included in each PCR assay. The positive control is a DNA template obtained from an *E. coli* strain possessing the target gene tested, while a negative control is the DNA from a non-pathogenic *E. coli* isolate (no virulence genes harbored, not belonging to the 5 serogroups) and the other is constituted by a sample without the template added.

The Real-Time PCR procedure requires an inhibition/extraction control. In particular, two different internal amplification controls (IACs) can alternatively be used:

- Commercially available IACs, often included in the kit with the mastermix for Real Time PCR.
- The open formula pUC 19 based internal amplification control IAC developed by Fricker *et al.* (2007. *Appl Environ Microbiol* 73:1892-1898). Approximately 100 copies of target DNA (pUC 19) should be used per PCR reaction. The size of the IAC was 119 bp. The latter may be also used as an extraction control by adding 100 copies of the pUC 19 plasmid to the sample aliquot, prior to the DNA purification step.

References

- Fricker, M., Messelhusser, U., Busch, U., Scherer, S., Ehling-Schulz, M. Diagnostic real-time PCR assays for the detection of emetic *Bacillus cereus* strains in foods and recent food-borne outbreaks. *Appl. Environ. Microbiol.* 2007, 73, pp. 1892–1898
- Perelle S, Dilasser F, Grout J, Fach P. Detection by 5'-nuclease PCR of Shiga-toxin producing *Escherichia coli* O26, O55, O91, O103, O111, O113, O145 and O157:H7, associated with the world's most frequent clinical cases. *Mol Cell Probes* 2004; 18:185-92.
- Perelle S, Dilasser F, Grout J, Fach P. Detection of *Escherichia coli* serogroup O103 by real-time polymerase chain reaction. *J Appl Microbiol.* 2005;98(5):1162-8.
- United States Department of Agriculture. Primer and Probe Sequences and Reagent Concentrations for non-O157 Shiga Toxin-Producing *Escherichia coli* (STEC) Real-Time PCR Assay. MLG 5B Appendix 1.01, https://www.fsis.usda.gov/wps/wcm/connect/0330211c-81ab-4e97-a9f3-d425f5759ee1/MLG_5B_Appendix_1_01.pdf?MOD=AJPERES