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Allegato 1

European Union Reference Laboratory for *E.coli*
Department of Veterinary Public Health and Food Safety
Unit of Foodborne Zoonoses and Veterinary Epidemiology
Istituto Superiore di Sanità

4th Inter-laboratory study on Verocytotoxin-producing *E. coli* (VTEC): detection and identification in milk samples

Laboratory procedure

The method to be used is the ISO-CEN draft Technical Specification “Horizontal method for the detection of Shiga toxin-producing *Escherichia coli* (STEC) belonging to O157, O111, O26, O103 and O145 serogroups - Qualitative Real-time polymerase chain reaction (PCR)-based Method”. This method has been developed by the “STEC *ad hoc* Group” within the WG 6 of Technical Committee 275 of the European Normalisation Committee (CEN TC275/WG6).

The method targets both virulence genes (*vtx1* and *vtx2*, and *eae*) and serogroup-specific genes for O26, O103, O111 and O145. It shall be used for screening samples, followed by the isolation of the VTEC strains responsible for the positive PCR reactions.

The method is applied to an enrichment culture performed by adding a test portion of 25 ml of the milk sample to 225 ml of mTSB containing 12 mg/l Acriflavin, and incubating for 18 -24 h at 37°C ± 1°C (see **Flow diagram** below). One ml aliquot of such a culture is used for DNA extraction and purification.

This step shall be accomplished according to the ISO 20837:2006 “Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of foodborne pathogens - Requirements for sample preparation for qualitative detection”.

The remaining culture shall be stored at 4°C for the isolation steps that will follow a positive PCR result.

Real-time PCR is performed using the primers and probes described in the included as **Annex 1**. Amplification conditions to be applied will depend on the system used, and will refer to the instructions supplied with the instrument and kit of choice.

The method is sequential:

Step 1: Enrichment of the sample.


Step 2: Detection of the genes *vtx1*, *vtx2* and *eae*.

Step 3: Samples positive for both *vtx* and *eae* at the second step are tested for the serogroup-associated genes (molecular serogrouping).

Step 4: Isolation of the VTEC strain; samples positive at the same time for *vtx*, *eae*, and serogroup-associated genes are submitted to a further step aimed at isolation of the VTEC strain. This requires serogroup-specific enrichment based on IMS or other immuno-capture suitable approaches. A guideline for the isolation of the different VTEC serogroups is included as **Annex 2**.

Step 5: Characterization of the isolate i.e. identification, detection of *vtx* genes, the *eae* gene and the serogroup gene.

Flow-diagram for the detection of VTEC non-O157 in milk samples

Annex 1

Real-time PCR for detection and identification of VTEC

1. Principle of the method

This Real-time PCR protocol aims at the detection and identification of the major virulence genes and serogroup-associated genes characterising the VTEC strains that are considered to be pathogenic to humans. These genes include:

- 1) *vtx* genes (*vtx1*, *vtx2* and its variants) encoding the Verocytotoxins (Shiga toxins), the main virulence factors of VTEC;
- 2) the *eae* gene, encoding a 90KDa protein, the intimin, which is the key factor for the induction of the “attaching and effacing” lesion on the enterocyte, a typical feature of the pathogenic VTEC strain;
- 3) genes associated with the VTEC serogroups that are mainly isolated from human cases of severe disease: O157, O26, O111, O103, and O145. The genes are either comprised in the operons encoding the different lipopolysaccharides (LPS) constituting the O antigens or are anyhow associated to each serogroup in a unique manner (Table 1).

2. Operating procedures

The protocol is based on the 5' nuclease PCR assay. Considering that Real-time PCR may use different instruments and probes labelling chemistry, the amplification conditions to be applied may vary depending on the system used. Refer to the instructions supplied with the instrument and kit of choice.

The primers and probes to be used are listed in the tables below. The chemistry of the reporter and quencher fluorophores are not indicated being largely dependent on the Real-time PCR systems available in each laboratory.

3. Controls

A VTEC strain harbouring *vtx1*, *vtx2* and *eae* genes should be used as positive control for the virulence genes. An example is the *E. coli* O157 EDL933 reference strain (ATCC no 43895). For serogroup-associated genes, reference strains belonging to each of the serogroup should be used.

The Real-time PCR procedure requires an **inhibition/extraction control**. Details on the possible systems to be used as inhibition/extraction control are given below. In particular, two different internal amplification controls (IACs) can alternatively be used:

- A commercially available TaqMan[®] Exogenous Internal Positive Control (Applied Biosystems, Foster City, CA, USA). The kit includes all reagents necessary (primers, a Vic[™] probe, IAC target DNA and blocking solution). The IAC target DNA must be diluted 10 times to achieve a copy number of approximately 100 per PCR reaction. The PCR product length is not declared to the customer.
- The open formula pUC 19 based internal amplification control IAC developed by Fricker et al. (2007. Appl Environ Microbiol 73, 1892-1898). Approximately 100 copies of target DNA (pUC 19) should be used per PCR reaction. The size of the IAC is 119 bp.

The latter may be also used as an extraction control by adding 100 copies of the pUC 19 plasmid to the sample aliquot prior to the DNA purification step.

Table 1. Degenerate primers and TaqMan probes used in 5' nuclease PCR assays for the detection of virulence genes

Target gene (Ref.)	Forward primer, reverse primer and probe sequences (5'-3') ^a	Amplicon size (bp)	Location within sequence	GenBank accession number
<i>vtx1</i> (1)	TTTGTYACTGTSACAGCWGAAGCYTTACG CCCCAGTTCARWGTRAGRTCMACRTC Probe- CTGGATGATCTCAGTGGGCGTTCTTATGTAA	131	878–906 983–1008 941–971	M16625
<i>vtx2</i> (1)	TTTGTYACTGTSACAGCWGAAGCYTTACG CCCCAGTTCARWGTRAGRTCMACRTC Probe- TCGTCAGGCACTGTCTGAAACTGCTCC	128	785–813 785–813 838–864	X07865
<i>Eae</i> (2)	CAT TGA TCA GGA TTT TTC TGG TGA TA CTC ATG CGG AAA TAG CCG TTA Probe- ATAGTCTCGCCAGTATTCGCCACCAATACC	102	899-924 1000-979 966-936	Z11541

^a In the sequence Y is (C, T), S is (C, G), W is (A, T), R is (A, G), M is (A, C)

(1) Perelle S. et al. Mol Cell Probes 2004 18: 185–192

(2) Møller Nielsen E. and Thorup Andersen M. J Clin. Microbiol. 2003 41: 2884-2893

(3) Perelle S. et al. J. Appl. Microbiol. 2005 98: 1162–1168

Table 2. Primers and probes used for the amplification of O antigen-specific genes in 5' nuclease PCR assays

Target gene (serogroup)	Forward primer, reverse primer and probe sequences (5'-3')	Amplicon size (bp) (Ref.)	Location within sequence	GenBank accession number
<i>§rfbE</i> (O157)	TTTCACACTTATTGGATGGTCTCAA CGATGAGTTTATCTGCAAGGTGAT Probe- AGGACCGCAGAGGAAAGAGAGGAATTAAGG	88 (1)	348–372 412–435 381–410	AF163329
<i>§wbdI</i> (O111)	CGAGGCAACACATTATATAGTGCTTT TTTTTGAATAGTTATGAACATCTTGTTTAGC Probe- TTGAATCTCCCAGATGATCAACATCGTGAA	146 (1)	3464– 3489 3579– 3609 3519– 3548	AF078736
<i>§wzx</i> (O26)	CGCGACGGCAGAGAAAATT AGCAGGCTTTTATATTCTCCAACCTTT Probe- CCCCGTAAATCAATACTATTTACGAGGTTGA	135 (1)	5648– 5666 5757– 5782 5692– 5724	AF529080
<i>§jhp1</i> (O145)	CGATAATATTTACCCCACCAGTACAG GCCGCCGCAATGCTT Probe- CCGCCATTCAGAATGCACACAATATCG	132 (1)	1383– 1408 1500– 1514 1472– 1498	AF531429
* <i>wzx</i> (O103)	CAAGGTGATTACGAAAATGCATGT GAAAAAAGCACCCCGTACTTAT Probe- CATAGCCTGTTGTTTTAT	99 (3)	4299– 4323 4397– 4375 4356– 4373	AY532664

(1) Perelle S. et al. Mol Cell Probes 2004 18: 185–192

(2) Møller Nielsen E. and Thorup Andersen M. J Clin. Microbiol. 2003 41: 2884-2893

(3) Perelle S. et al. J. Appl. Microbiol. 2005 98: 1162–1168

Annex 2:

Isolation of VTEC O26, O103, O111 AND O145 from Real-time PCR positive samples

The Real-time PCR method proposed for the detection of VTEC belonging to the pathogenic serogroups O26, O103, O111 and O145 is based on screening samples for the presence of *vtx* (VT-coding genes), *eae* (intimin-coding gene) and specific serogroup-associated genes. Samples positive at the same time for *vtx*, *eae*, and a serogroup-associated gene must be submitted to a further step aimed at the isolation of the VTEC strain, to confirm the simultaneous presence of the genes in the same live bacterial cell.

Flow of the operations include:

- 1) Submit the enrichment culture maintained as a back-up to serogroup-specific enrichment (SSE), using immuno-capture systems such as immuno-magnetic separation (IMS) or equivalent following instructions supplied by the manufacturer.
- 2) Streak the SSE onto Tryptone-bile-glucuronic medium (TBX) or onto a specific selective medium where available (see Note 1) and incubate for 18 to 24 hours at 37°C.
- 3) Pick from 10 to 50 colonies with *E. coli* morphology or with characteristic aspect according to the medium used (see Note 1) and point-inoculate on nutrient agar (NA) (see Note 2) and H₂O (the colonies may be pooled in water up to a number of ten per pool).
- 4) Perform conventional PCR or Real-Time PCR on the H₂O pools to assess the presence of the *vtx* and the *eae* genes.
- 5) If a PCR pool is positive, go back to NA and assay the individual colonies forming the positive pool in order to select one single positive colony.
- 6) Identify the isolates as *E. coli* and confirm the serogroup the sample was positive for in the screening PCR assay (see Note 3).

NOTE 1: For VTEC O26 isolation, a differential solid media (MacConkey) containing Rhamnose instead of lactose is commercially available (RMAC). It is very effective in distinguishing VTEC O26 strains, which do not ferment Rhamnose, from other *E. coli*. Differential and confirmation plating media for VTEC serotypes O26, O103, O111, O145, and sorbitol-positive and -negative O157 have been recently proposed (1, 2). Enterohaemolysin Agar can also be used. It detects Enterohaemolysin production, which is a common feature of VTEC pathogenic to humans (3).

NOTE 2: There are several types of nutrient agar media available commercially either ready-to-use plates or prepared in house from dehydrated powders. Every type of non-selective nutrient agar media (e.g. TSA), including Enterohaemolysin Agar, is suitable for the purpose of maintaining the colonies for further characterisation.

NOTE 3: Colony confirmation as *E. coli* may be achieved by using commercial biochemical galleries or by assessing the indole production. Confirmation of the serogroup the sample was positive to in the screening PCR assay may be achieved either by PCR or by agglutination with commercial antisera.

References

Posse B, De Zutter L, Heyndrickx M, Herman L, 2008b. Quantitative isolation efficiency of O26, O103, O111, O145 and O157 STEC serotypes from artificially contaminated food and cattle faeces samples using a new isolation protocol. *J Appl Microbiol*, 227-235.

Possé B, De Zutter L, Heyndrickx M, Herman L, 2008. Novel differential and confirmation plating media for Shiga toxin-producing *Escherichia coli* serotypes O26, O103, O111, O145 and sorbitol-positive and -negative O157. *FEMS Microbiol Lett* 282,124-131.

Beutin L, Zimmermann S, Gleier K. Rapid detection and isolation of shiga-like toxin-producing *Escherichia coli* by direct testing of individual enterohemolytic colonies from washed sheep blood agar plates in the VTEC-RPLA assay. *J Clin Microbiol*. 1996; 34:2812-4.

Figure 1. Flow-diagram of the isolation procedure


Figure 1. Flow-diagram of the isolation procedure