



5th annual workshop of the national reference laboratories for *E. coli* in the EU



STEC surveillance in France:
official laboratory network
and results of 2010 ring trial



Dr. Estelle LOUKIADIS*

8th October 2010, ISS, Rome, Italy

VetAgro Sup

* *VetAgro Sup Campus Vétérinaire de Lyon, UMAP*

e.loukiadis@vetagro-sup.fr



STEC surveillance in France: official laboratory network and results of 2010 ring trial



Plan

1- Official laboratory network


- Official controls organisation
- Networks creation/extension in 2009 and 2010
- Territorial distribution

2- Results of 2010 ring trial

Dr. Estelle LOUKIADIS

5th annual workshop of the NRL for *E. coli* in the EU, 8th October 2010





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Official controls organisation

→ Annual studies

(Directive EU 2003/99/CE; Regulation EU N° 882/2004)

Objective: Estimate STEC* prevalence in risk-associated food
(* *E. coli* possessing *eae* and *stx* genes and belonging to the 5 major EHEC-associated serogroupes *i.e.* O157, O26, O103, O111 and O145)



Dr. Estelle LOUKIADIS

5th annual workshop of the NRL for *E. coli* in the EU, 8th October 2010



STEC surveillance in France: official laboratory network and results of 2010 ring trial

Year	Matrix	Nb of samples screened	Nb of strains isolated	STEC * prevalence
2007	Ground beef batch (+4°C)	3605	11 STEC* (5 <i>aeae+</i> <i>E. coli</i>)	11/3605 0,3% CI ₉₅ [0,2 – 0,5%]
	Raw milk cheese	392	0	0/392 ≤ 0,9%
2008	Ground beef trimming (-20°C)	992 (4 units/samples)	10 STEC* (13 <i>aeae+</i> <i>E. coli</i>)	10/992 1% CI ₉₅ [0,6 – 1,8%]
2009	Minced beef at retail	1527	2 STEC* (3 <i>aeae+</i> <i>E. coli</i>)	2/1527 0,1% CI ₉₅ [0,04 – 0,5%]
	Raw milk cheese	1911	17 STEC* (14 <i>aeae+</i> <i>E. coli</i>)	17/1911 0,9% CI ₉₅ [0,6 – 1,4%]
2010	Minced beef at retail	1403 out of 2500	4 STEC* (9 <i>aeae+</i> <i>E. coli</i>)	3/1403



STEC surveillance in France: official laboratory network and results of 2010 ring trial



Official controls organisation

Annual studies


→ Samples analyses

- Official method: CEN/ISO TS





« Microbiology of food and animal feeding stuffs - Horizontal method for the detection of STEC belonging to O157, O111, O26, O103 and O145 serogroups – Qualitative method »





STEC surveillance in France: official laboratory network and results of 2010 ring trial





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Official controls organisation


Annual studies

→ Samples analyses

- First screening by official laboratories
- Confirmation by the NRL




STEC surveillance in France: official laboratory network and results of 2010 ring trial





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Plan

1- Official laboratory network

- Official control organisation
- Networks creation/extension in 2009 and 2010
- Territorial distribution

2- Results of 2010 ring trial

Dr. Estelle LOUKIADIS

5th annual workshop of the NRL for *E. coli* in the EU, 8th October 2010



STEC surveillance in France: official laboratory network and results of 2010 ring trial

Network creation/extension in 2009 and 2010



French regulation

(Regulation EU N° 882/2004)

- Publication of a call for application
- General requirements
(accreditation/information exchange)
- Specific requirements
(technical accreditation/material/training/ring trial)




➔ 2 networks created (O157 and STEC)

Dr. Estelle LOUKIADIS

5th annual workshop of the NRL for *E. coli* in the EU, 8th October 2010





STEC surveillance in France: official laboratory network and results of 2010 ring trial




Network creation/extension in 2009 and 2010



Role of the NRL

- Provide an official method
- Provide technical support for the call and for the analyses of the laboratories' files
- Organize trainings
- Organize a ring trial





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Plan

1- Official laboratory network

- Official control organisation
- Networks creation/extension in 2009 and 2010
- Territorial distribution

2- Results of 2010 ring trial

Dr. Estelle LOUKIADIS


5th annual workshop of the NRL for *E. coli* in the EU, 8th October 2010




STEC surveillance in France: official laboratory network and results of 2010 ring trial

Territorial distribution

➔ **O157 network:**
23 laboratories (●)





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Territorial distribution

➔ **STEC network:
17 laboratories (•)**

...Probably
19 laboratories in 2011





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Plan

1- Official laboratory network

- Official control organisation
- Networks creation/extension in 2009 and 2010
- Territorial distribution

2- Results of 2010 ring trial

Dr. Estelle LOUKIADIS

5th annual workshop of the NRL for *E. coli* in the EU, 8th October 2010





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Results of 2010 ring trial



Objectives:


- Quality assurance
- Extension of the STEC network in 2010

(Participation is mandatory in the call for application of French ministry of agriculture)





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Results of 2010 ring trial



Method:

- Target: Detection of *eae* and *stx* genes in minced beef

- Official method:

Detection of *stx* and *eae* genes by RT-PCR with internal control

(1 method provided by the NRL adapted from CEN/ISO TS)



Dr. Estelle LOUKIADIS

5th annual workshop of the NRL for *E. coli* in the EU, 8th October 2010





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Results of 2010 ring trial

→ Participants:

- Food inspections laboratories in France
- Only candidates that were first selected by the French ministry of Agriculture (general and specific requirements indicated in the call)
- 19 laboratories participate + the NRL



STEC surveillance in France: official laboratory network and results of 2010 ring trial

Results of 2010 ring trial


→ Sample preparation:

- Minced meat (10 kg)
(collected at retail)
- Enrichment performed by the NRL (37°C – 20h)
(enrichment broth checked for the absence of *stx* and *eae* genes)
- Artificial contamination (or not) of enrichment broth
(level of contamination checked before and after contamination)
- 8 samples (5ml) of enrichment broth
(same test portion for each participant)



Dr. Estelle LOUKIADIS

5th annual workshop of the NRL for *E. coli* in the EU, 8th October 2010



STEC surveillance in France: official laboratory network and results of 2010 ring trial




Results of 2010 ring trial

➔ Characteristics of the 8 samples sent


Strains /samples	A	B	C	D	E	F	G	H
<i>E. coli</i> O157:[H7] <i>stx</i> ₁ ⁺ , <i>stx</i> ₂ ⁺ , <i>eae</i> ⁺	5x10 ⁴ CFU/ml						5x10 ⁴ CFU/ml	
<i>E. coli</i> O157:[H7] <i>stx</i> ₁ ⁻ , <i>stx</i> ₂ ⁺ , <i>eae</i> ⁺		10 ⁵ CFU/ml						
<i>E. coli</i> O26:[H11] <i>stx</i> ₁ ⁻ , <i>stx</i> ₂ ⁻ , <i>eae</i> ⁺				10 ⁵ CFU/ml				
<i>E. coli</i> O26:[H11] <i>stx</i> ₁ ⁺ , <i>stx</i> ₂ ⁻ , <i>eae</i> ⁺					10 ⁵ CFU/ml			
<i>E. coli</i> O139:H – <i>stx</i> ₁ ⁻ <i>stx</i> ₂ ⁺ , <i>eae</i> ⁻			10 ⁵ CFU/ml					

stx⁺ or *eae*⁺ samples:
Detectability

stx⁻ or *eae*⁻ samples:
Specificity



STEC surveillance in France: official laboratory network and results of 2010 ring trial



Results of 2010 ring trial

➔ Detection of *stx* genes by RT-PCR

Détection des gènes	Echantillons	Valeur attendue	Laboratoires participants																			
			L1	L2	L3	L4	L5	L6	L7	L8	L9	L10	L11	L12	L13	L14	L15	L16	L17	L18	L19	L20
<i>stx</i>	A	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
	B	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
	C	+	+	+	+	+	+	+	+	+	+	+	-	+	+	+	+	+	+	+	+	+
	D	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	E	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
	F	-	-	-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	G	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
	H	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

Detectability (samples A, B, C, E, G)

19 laboratories do detect all the *stx*+ samples

1 laboratory fails to detect *stx* gene in sample C (internal control OK)

Specificity (samples D, F, H)

19 laboratories get expected results (*stx*-)

1 laboratory gets false positive results (cross contamination?)

STEC surveillance in France: official laboratory network and results of 2010 ring trial

Results of 2010 ring trial

➔ Detection of *eae* genes by RT-PCR

Détection des gènes	Echantillons	Valeur attendue	Laboratoires participants																			
			L1	L2	L3	L4	L5	L6	L7	L8	L9	L10	L11	L12	L13	L14	L15	L16	L17	L18	L19	L20
<i>eae</i>	A	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
	B	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
	C	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
	D	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
	E	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
	F	-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-	-	-	-
	G	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
	H	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-


Detectability (samples A, B, D, E, G)

All laboratories (20) do detect all the *eae* + samples

Specificity (samples C, F, H)

19 laboratories get expected results (*eae*-)

1 laboratory gets false positive results (cross contamination ?)



STEC surveillance in France: official laboratory network and results of 2010 ring trial



Results of 2010 ring trial



Conclusion

- Number of samples screened and statistical interpretation!

Ex: Detectability (evaluation with 5 samples)

If 5 expected results/5, discrepancy frequency CI_{95} [0,48 – 1]

If 4 expected results/5, discrepancy frequency CI_{95} [0,28 – 0,99]

Lamy *et al.*, Simple table for estimating confidence interval of discrepancy frequencies in microbiological safety evaluation. Journal of microbiology methods (2004) 137-139.




Results of 2010 ring trial

➔ Conclusion

- 18 out of 20 laboratories get all the expected results
- 2 laboratories failed

➔ New ring trial next week!





STEC surveillance in France: official laboratory network and results of 2010 ring trial



Delphine Thevenot



Christine
Mazuy



Christine
Bavai



Lysiane
Dunière



Fanny
Savoye



Stéphane
Mischczycha



Sarah
Ganet



Marion
Bouvier




Audrey
Gleizal



Franck
Ferre

And Manuelle Neto, Cécile Dutrieux, Françoise Clavel, Patrice Zurcher!!



STEC surveillance in France: official laboratory network and results of 2010 ring trial



Thank you
for your attention!



VetAgro Sup

VetAgro Sup Campus vétérinaire de Lyon
UMAP

1 Avenue Bourgelat
69380 Marcy l'Etoile FRANCE

Tel : +33 4 78 87 25 53 Fax: +33 4 78 87 26 69